

Sizes

The Qseals QFMS06 is available in at least 9 different sizes:

D1	Unit
25.00	mm
30.00	mm
33.00	mm
35.00	mm
38.00	mm
45.00	mm
65.00	mm
70.00	mm
90.00	mm

Technical specifications

Temperature	... till ... °C
Pressure	... bar
Rotational velocity	... m/s
Rotate left	No
Rotate right	No
Balanced	No

Alternatives

The Qseals QFMS06 is an alternative to at least 52 products of other brands:

- AES M03U
- Anga VT EN 12756
- Burgmann BT-C56.KU
- Flowserve M412LZ199
- Flowserve M412LZ426
- Flowserve M412LZ416
- Flowserve M412LZ414
- Flowserve M412LZ389
- Flowserve M412LZ313
- Flowserve M412LZ303
- Flowserve M412LZ270
- Flowserve M412LZ268
- Flowserve M412LZ219
- Flowserve M412LZ191
- Flowserve M412LZ192
- Flowserve M412LZ173
- Flowserve M412LZ125
- Flowserve M412LZ124
- Flowserve M412LZ112
- Flowserve M412LZ061
- Flowserve M412LZ043
- Flowserve M412LZ042
- Flowserve M412LZ019
- Flowserve M411LZ733
- Flowserve M412LZ437
- And 27 others...

Applications

The Qseals QFMS06 is used in at least 6 pumps and other applications:

- Mitsubishi CVD-25L

- Viking LVP19
- Viking LVP23
- Viking LVPO5
- Viking LVPO8
- Viking SG414