

Sizes

The Qseals QFMS19 is available in at least 3 different sizes:

D1	Unit
24.00	mm
32.00	mm
38.00	mm

Technical specifications

Temperature	... till ... °C
Pressure	... bar
Rotational velocity	... m/s
Rotate left	No
Rotate right	No
Balanced	No

Alternatives

The Qseals QFMS19 is an alternative to at least 2 products of other brands:

- AES M010SL
- Flowserve Europac 306

Applications

The Qseals QFMS19 is used in at least 17 pumps and other applications:

- Nowa 4016BN138422
- Sterling SIHI CEHA 6101
- Sterling SIHI CEHA 6107
- Sterling SIHI CEHA 6106
- Sterling SIHI CEHA 6105
- Sterling SIHI CEHA 6104
- Sterling SIHI CEHA 6103
- Sterling SIHI CEHA 6102
- Sterling SIHI CEHA 5108
- Sterling SIHI CEHA 5101
- Sterling SIHI CEHA 5107
- Sterling SIHI CEHA 5106
- Sterling SIHI CEHA 5105
- Sterling SIHI CEHA 5104
- Sterling SIHI CEHA 5103
- Sterling SIHI CEHA 5102
- Sterling SIHI CEHA 6108